luminette® Glasses 3

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Model: Glasses 3
EN - ES - FR - DE - NL

User Manual Luminette 3
EU-US V1


Figure 1 – The Luminette and its accessories

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



Figure 3 – Positioning while wearing eyeglasses

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
EN

I. Introduction

Congratulations! You've just bought a Luminette®. Welcome to the world of portable light therapy.

Before using your Luminette®, please read this notice carefully and keep it somewhere safe. Also, please take the time to familiarise yourself with your new appliance using Figure 1.

A. The role of light

The human brain uses the light/dark cycle to synchronise internal functions including mood, appetite, digestion, sleep, libido, etc. These circadian rhythms oscillate over a period of 24 hours, forming your biological clock.

Melatonin is what is often called the sleep hormone. The brain starts producing it when the sun sets, in order to encourage us to sleep. But as the sun rises and the light increases, the brain stops producing it. We get up easily, full of energy.

As well as regulating our body clock, light also greatly affects mood, increases vigilance and improves physical and intellectual performance.

B. Light therapy

In certain situations, our body clock is no longer in phase with the external environment. The circadian rhythms are out of sync.

This is particularly true in winter when there is insufficient light; some people develop winter blues. This can also happen in people who go to sleep too early or too late. Jet lag is another form of a disrupted body clock. Finally, problems relating to shift work are now recognised as real health problems.

All these problems can be relieved with the correct use of light and dark. When light is not naturally available, light therapy is the solution which offers the most advantages for the least unwanted side effects.

C. Applications

- Winter blues
- Fatigue - lack of energy
- Circadian rhythm sleep disorder

- Jet lag
- Working antisocial hours

The Luminette is intended to make people feel more energetic, to improve sleeping patterns, to increase alertness, to improve intellectual performance, and to provide relief from the winter blues. If you suffer from depression, consult a professional and do not use the Luminette for self-treatment.

This product is not a medical device.

D. Contraindications

No matter your reason for using Luminette®, we recommend not exceeding the recommended daily duration of use.

There are very few contraindications.

If you have an eye disease (glaucoma, retinitis or retinopathy, or macular degeneration), we recommend that you consult a specialist before using the Luminette®. If you have cataracts, light therapy might be less effective.

For people with a bipolar disorder, the Luminette® could be beneficial if used at the right time when someone is in a depressive phase. But care should be taken as there is always the risk of triggering a manic phase. It is therefore recommended to time a session for midday, when the person is least sensitive to light and therefore less likely to swing into a manic phase. The same is true for cases of hysteria, neurosis and psychosis.

It is also not advisable to use the Luminette® whilst driving a vehicle, or using machinery requiring concentration and a wide field of vision.

E. Possible reactions

Possible adverse reactions to light from this appliance are relatively mild and of a transient nature. They may include:

- Headache.
- Eyestrain or visual disturbances, intraocular pressure.
- Nausea.
- Agitation, sedation, jumpiness/jitteriness and very rarely hypomania or mania.

Measures to avoid, minimise or alleviate these adverse reactions:

- Always use the appliance in a well-lit room.
- Reduce the light intensity if these reactions occur.
- Stop using the appliance for a few days to make the reactions disappear and then try again.

F. Safety

This appliance conforms to technical regulations and safety requirements applicable to electrical equipment. All repairs and maintenance must be carried out by a specialist. Ask your supplier for advice.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

To ensure the longevity of your Luminette®:

- Only use original accessories.
- Never immerse your Luminette® in water.
- Don't wash your Luminette® in running water.
- Don't leave your Luminette® in the sun, or in damp or dusty conditions.

Attention!

- No modification to the device or its accessories is authorised.
- This device cannot be attached to an extension cord or a multi-socket adapter.
- Always use the device with the adaptor provided or any other compatible system as specified in the manual.
- The adaptor constitutes the network dividing device and must always be easily accessible.
- Be aware that the external envelope of the central case can reach a temperature of 49°C.
Applied parts temperatures exceeds 41°C which are safe for use for the maximum duration of a session (i.e. 45 min)
- In case of mutual interference, electromagnetic or other with other devices, move the device away.

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by Lucimed could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Luminette3, including cables specified by Lucimed. Otherwise, degradation of the performance of this equipment could result.

II. Use

A. When can I use my Luminette®?

1. To fight the winter blues

You lack energy and motivation, you have difficulty getting up in the morning, you feel like hibernating, you snack more, your libido is falling? You could be suffering from winter blues, which results from lower light levels in autumn and winter.

We advise a light therapy session as soon as you wake up in the morning.

2. To fight fatigue and get energy back

Luminette® reproduces the beneficial power of daylight and stimulates special receptors in the eye to trigger the body's energising response to light. This helps you to fight energy dips, to feel more energetic, active and alert, and to improve your mood.

We advise a light therapy session as soon as you wake up.

3. To treat circadian rhythms sleep disorders

Some sleep problems result from a phase advance or delay in your biological clock.

In which case, light therapy could be the solution.

Your body clock is slow

Do you tend to fall asleep too late and have a hard time getting up in the morning?

We advise a light therapy session as soon as you wake up.

Your body clock is fast

Do you have a hard time staying awake in the evening and do you wake up too early in the morning? We advise a light therapy session at the end of the afternoon or the beginning of the evening, two or three hours before the time you plan to go to bed.

4. To treat jet lag

Jet lag is the result of crossing several time zones very rapidly, usually in a plane, and results in disturbed body rhythms: insomnia, drowsiness, physical and intellectual fatigue, mood swings, difficulty concentrating...

Some simple steps and the use of Luminette® can considerably reduce the symptoms.

- During the flight

Set your watch to the time zone of your destination, drink lots of water and try not to drink alcohol, coffee or tea.

- When you arrive

If you are travelling east, when you arrive we advise you to avoid sunlight on the first morning (by wearing blue blocker glasses or high sun protection sunglasses) and to get as much as you can in the afternoon (by, for example, using your Luminette®) in order to maintain low levels of melatonin. During the following days, put your session forward by two hours a day, taking care not to go past your normal bedtime.

If you are travelling west, once you arrive, we advise you to stay awake as long as it is still light, and to wear your Luminette® at the beginning of the afternoon (between midday and 3pm). But as soon as it is dark, go to sleep.

5. To optimise shift work

If you work for more than three nights in a row and you have a hard time staying awake at work and/or on your way home, then we recommend a light therapy session just before starting work and, if you get drowsy, a second session during your work break.

We also advise that you talk to your GP.

Would you like a personalised programme? Download the free MyLuminette app available on App Store and Google Play Store.

B. How long does a light therapy session last?

The recommended duration of a session depends on the intensity of the light received by the eye.

Your Luminette® allows you to choose from three different light intensities.

The table below will tell you the recommended duration of use for your Luminette® depending on the position you choose.

Position (when turning on appliance)	Recommended duration of session (in minutes)
1	45
2	30
3	20

The chosen light intensity is indicated by the 3 light intensity indicators (see figure 1). Once the session time is up, the glasses will blink three times and automatically turn off.

C. How often do I have to use my Luminette?

Use the Luminette every day at the right time to achieve the desired results. Except when you try to beat jet-lag, it is advised to use the Luminette at the same time every day.

D. When do I start to notice improvement ?

Some users start to feel the energising effect after the first use. Others feel it after 1 or 2 weeks of daily use.

E. How to charge my Luminette ?

Before using your Luminette® for the first time, we recommend to charge it.

To do so:

- Connect the end of the micro-USB cable to the slot in your Luminette® (see figure 1) and plug the other end (a standard USB) in the charger.
- Plug the charger into the mains.
- When the device is completely charged (the charge time is 2 hours), unplug the charger.

You can also charge your Luminette® by connecting it to your computer using the USB cable supplied.

Charge indicator:

- When your Luminette® is charging, the 3 battery indicators will continually flash slowly in a cycle.
- When your Luminette® is fully charged, these 3 battery indicators stop flashing and remain lit.
- If there is an error (battery or circuit failure), one of the indicators flashes rapidly without interruption.

The device may not be used when it is charging.

Battery life:

The battery life will depend on the intensity you choose. A correctly charged Luminette® device should lit for at least 10 sessions.

When you turn the device off, between 1 and 3 battery indicators light up to give you an indication of the remaining battery life (3 indicators = the battery is fully charged; 2 indicators = the battery is 2/3 charged; 1 indicator = the battery is 1/3 charged; 1 flashing indicator = the battery is insufficiently charged to be able to complete a session).

We recommend recharging your Luminette® when there is only one indicator remaining when you turn the device off.

F. Switching on and off

Attention: before your first use, be sure you take off the protective film over the

hologram.

Ensure that the battery is charged so that your Luminette® will work for the entire session chosen. Press the Start-Stop button briefly to turn your Luminette® on (see figure 1).

To change the intensity of the light, just press the Start-Stop button briefly. The intensity will change in a cycle: level 1, then level 2, level 3, and back to level 1, etc.

Luminette remains on for a set period, blink three times when the period is over and automatically shut down. To turn your Luminette® off in the middle of a session, hold down the On-Off button (for 2 seconds).

G. How to wear my Luminette®

For the best results, you need to wear your Luminette® correctly.

To properly position your Luminette®:

- Set the nosepiece into one of the three slots (see Figure 1). A “click” sound means that your nosepiece is correctly inserted. To take it out, pull firmly on it.
- Place the padded part of the nosepiece where your nose and forehead meet. The lower edge of the holographic visor should be located just above your line of sight (see figure 2) when you look into a mirror. Make sure the blue light strip sits nicely across your eyes.

Please note that if you wear glasses, the nose bracket of the Luminette® should be behind your glasses (see figure 3).

Once the session has begun, keep your eyes open. Carry on with whatever you are doing. Do not look directly into the LEDs.

Find the perfect fitment : Luminette is one-size-fits-all; however, if it needs tweaking for a tighter or looser fit, bend gently the temples in or out. For an even more secure fit, insert the supplied white plastic “ear-hooks” at the end of the temples.

H. How to look after my Luminette®

This device is meant exclusively for personal use. Between each use, store the device in its protective case in a dry place.

To clean your Luminette®, use a soft cloth moistened with a little clean water. Be careful to protect the fragile holographic zone; avoid touching it with your

fingers. If, despite everything, it gets dirty, clean with cotton wool.

III. Myluminette app - instructions

A. Installation of the app

Search for the app called "Myluminette" in the Apple store or in Google Play store. Download and install it on your smartphone.

B. Create your profile

Once you have installed the app "Myluminette", you are asked to create a profile and enter some basic personal information.

C. Pair the app with my Luminette

Press and hold the on/off button of your Luminette for more than 5 seconds to start the pairing mode. The 3 LED light indicators will start flashing slowly when the pairing mode is activated.

Turn on the Bluetooth functionality on your smartphone and select Luminette in the list of detected Bluetooth-emitting devices. Your Luminette should be now paired with your smartphone. You can use the app and have access to your personal usage history.

IV. Guarantee and After Sales Service

A. Guarantee

Lucimed guarantees the operation of its appliances when used in accordance with the Instructions for Use, against all manufacturing defects. This guarantee covers spare parts and labour. If the appliance ceases to function, Lucimed guarantees free repair for 2 years from the date of purchase.

The guarantee is only valid when the sales receipt or the completed and date-stamped guarantee card is provided. The repair must be carried out within the guarantee period.

The guarantee does not cover any faults relating to accidental damage, incorrect use or negligence, or to modifications or repairs carried out by non-qualified personnel.

In no way does this guarantee give the user rights to damages.

Lucimed provides after sales service during and post the guarantee period in all countries where the appliance is officially sold and distributed.

In case of any problem, please visit www.myluminette.com for assistance.

Include the serial number and date of purchase on all correspondence.

Please note: if you ever lose these instructions for use, don't forget that you can download a new copy from www.myluminette.com.

B. Accessories

The following accessories are available upon purchase:

- Charging adaptor (EU): reference KYT050100BC-05
- Protective case: reference CASE02

V. Environment

To dispose of your appliance, do not include it with household waste but try to use a recycling centre in order to help protect the environment.

The rechargeable battery inside your Luminette® contains substances which could harm the environment.

So please remove the used battery from the appliance before disposing of it, and dispose of it separately in a battery recycling bin. If you can't remove the battery, please return the appliance to the vendor.

To remove the rechargeable battery:

- Ensure that the appliance is unplugged and that the battery is flat.
- Unscrew the two arms using a correctly sized screwdriver.
- Remove the plastic cover from the rechargeable battery.
- Unplug the rechargeable battery.

VI. Technical Specifications

A. Features of the appliance

General	
Trademark	Luminette®
Device model reference	Glasses 3
Adaptor model reference (EU) + cable	KYT050100BC-05
Electricity	
Nominal input voltage of the adapter	100-240V AC
Nominal input frequency of the adapter	50-60 Hz
Nominal output voltage of the adapter	5V DC
Nominal output power of the adapter	5W
Power type	
2 Rechargeable Lithium Ion battery	3.7V 150 mAh
Physical characteristics	
Dimensions	17 x 16 x 3 cm (6.7 x 6.3 x 1.2 inch.)
Weight	51 gr. (1.8 Oz)
Light intensity at Position 1	+ - 500 lux
Light intensity at Position 2	+ - 1,000 lux
Light intensity at Position 3	+ - 1,500 lux
Spectral characteristics	
Peak 1 Wavelength Bandwidth	468 nm 70 nm
Peak 1 Wavelength Bandwidth	570 nm 100 nm

Nominal Spectral irradiance (100% intensity)	21,8 $\mu\text{W}/\text{cm}^2$
Maximal Spectral irradiance (100% intensity)	25,4 $\mu\text{W}/\text{cm}^2$
Maximal variation of the output	$\pm 7\%$
Operating conditions	
Temperature	5 °C - +40 °C
Relative humidity	15 % - 93 %
Atmospheric pressure range	700 hPa - 1060 hPa
Storage conditions	
Temperature	-25 °C - +70 °C
Relative humidity	15 % - 93 %
Expected service life	
Device lifespan	6 years
Adaptor lifespan	3 years

B. Symbols

The following symbols may appear on the appliance:

Symbols	Description
	Manufactured for : Lucimed SA, Rue le marais 12a, 4530 Villers-le-Bouillet, Belgium Tel : +32 4 369 48 36
	Year of manufacture
	Consult the user manual supplied
	This device should not be disposed of in the normal waste stream but in an appropriate collection facility for recovery and recycling

	Carefully consult the user manual supplied
	The packaging of this device can be recycled
	CE stands for European Declaration of Conformity.

C. Manufacturer's liability

As manufacturers, we are only responsible for the safety, reliability and functioning of this appliance when it is used in accordance with the directions for use and when modifications and repairs are exclusively carried out by personnel approved by us.

Please contact the manufacturer if you require assistance with installing, using or maintaining your device, or if you encounter an unexpected problem.

Any serious incident involving Luminette must be notified to Lucimed and your country's qualified authority.

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